

## Product datasheet for **HP229000**

### **KIAA0701 (UHRF1BP1L) Human qPCR Primer Pair (NM\_015054)**

#### **Product data:**

Product Type:	qPCR Primer Pairs
Gene ID:	23074
Forward Sequence:	GCTTCACTCTCTGAAAACCTGGA
Reverse Sequence:	GAAACACAGAGAGGAGCCAGGT
ACCN:	<u><a href="#">NM_015054</a></u> , <u><a href="#">NM_015054.1</a></u> , <u><a href="#">BC127017</a></u> , <u><a href="#">BC003391</a></u> , <u><a href="#">BC014891</a></u> , <u><a href="#">BC026286</a></u> , <u><a href="#">BC026304</a></u> , <u><a href="#">BC062763</a></u> , <u><a href="#">BX647380</a></u> , <u><a href="#">NM_015054.2</a></u>
UniProt ID:	<u><a href="#">A0JNW5</a></u>
Synonyms:	SHIP164
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)