

Product datasheet for **HP229009**

CABP (CABP1) Human qPCR Primer Pair (NM_031205)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9478
Forward Sequence:	GAACTGCGAGATGCTTTCCGAG
Reverse Sequence:	CCATTGAGGTCCACATCTCGGA
ACCN:	<u>BC015006</u> , <u>NM_031205</u> , <u>NM_031205.1</u> , <u>NM_031205.2</u> , <u>NM_031205.3</u> , <u>BC030201</u> , <u>BC030201.1</u> , <u>BM716422</u> , <u>BQ188859</u> , <u>NM_031205.4</u>
UniProt ID:	<u>Q9NZU7</u>
Synonyms:	CALBRAIN; HCALB_BR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)