

Product datasheet for **HP229045**

TFII I (GTF2I) Human qPCR Primer Pair (NM_032999)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2969
Forward Sequence:	CGACTTTTGGCATTCCGAGGCT
Reverse Sequence:	GTGCTCTCCTTAATCGCCGTCT
ACCN:	<u>NM_032999</u> , <u>NM_032999.1</u> , <u>NM_032999.2</u> , <u>NM_032999.3</u> , <u>BC070484</u> , <u>BC070484.1</u> , <u>BC099907</u> , <u>BC099907.1</u> , <u>BC004472</u> , <u>BT007450</u>
UniProt ID:	<u>P78347</u>
Synonyms:	BAP135; BTKAP1; DIWS; GTFII-I; IB291; SPIN; TFII-I; WBS; WBSCR6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)