

Product datasheet for **HP229065**

HNRPDL (HNRNPDL) Human qPCR Primer Pair (NM_031372)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9987
Forward Sequence:	CAATAGTGCCTATGGTGGTGATC
Reverse Sequence:	GCCTTGCCATAAGTGCTCTGTTG
ACCN:	<u>NM_031372</u> , <u>NM_031372.1</u> , <u>NM_031372.2</u> , <u>NM_031372.3</u> , <u>BC011714</u> , <u>BC011714.2</u> , <u>BC007392</u> , <u>BC018961</u> , <u>BC071944</u>
UniProt ID:	<u>O14979</u>
Synonyms:	HNRNP; HNRPDL; JKTBP; JKTBP2; laAUF1; LGMD1G; LGMDD3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)