

Product datasheet for **HP229067**

PCDH9 Human qPCR Primer Pair (NM_203487)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5101
Forward Sequence:	CACTCCTCAGAAATGTCCCAGC
Reverse Sequence:	TGGTCTCCTAGACCACTGTCCAC
ACCN:	NM_203487 , NM_203487.1 , NM_203487.2 , BC047557 , BC136626 , BC136627 , BC144321 , BC150296 , BX537422
UniProt ID:	Q9HC56
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)