

## Product datasheet for **HP229160**

### FAM111A Human qPCR Primer Pair (NM\_198847)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63901
Forward Sequence:	GCACTGTGCAAGGATGGCAGAT
Reverse Sequence:	CCTCAACCTGAAAGTATCTGCCT
ACCN:	<u><a href="#">NM_198847</a></u> , <u><a href="#">NM_198847.1</a></u> , <u><a href="#">NM_198847.2</a></u> , <u><a href="#">BC071759</a></u> , <u><a href="#">BC071759.1</a></u> , <u><a href="#">BC013137</a></u> , <u><a href="#">BC041693</a></u> , <u><a href="#">BC054515</a></u> , <u><a href="#">BQ420366</a></u> , <u><a href="#">BQ948815</a></u> , <u><a href="#">BU428612</a></u> , <u><a href="#">NM_198847.3</a></u>
UniProt ID:	<u><a href="#">Q96PZ2</a></u>
Synonyms:	GCLEB; KCS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)