

## Product datasheet for **HP229164**

### RRAS2 Human qPCR Primer Pair (NM\_012250)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22800
Forward Sequence:	GAAGGACAACAGTTAGCACGGC
Reverse Sequence:	TCCGTGTTGGTTCTGGTGAAGG
ACCN:	<a href="#">NM_012250</a> , <a href="#">NM_012250.1</a> , <a href="#">NM_012250.2</a> , <a href="#">NM_012250.3</a> , <a href="#">NM_012250.5</a> , <a href="#">BC013106</a> , <a href="#">BC013106.1</a> , <a href="#">BG428983</a> , <a href="#">BP206021</a> , <a href="#">BP206589</a> , <a href="#">BQ020351</a> , <a href="#">BQ023701</a> , <a href="#">BX324335</a> , <a href="#">NM_012250.6</a>
UniProt ID:	<a href="#">P62070</a>
Synonyms:	NS12; TC21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)