

## Product datasheet for **HP229187**

### DBNDD2 Human qPCR Primer Pair (NM\_001048225)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55861
Forward Sequence:	CCTTGCGAAGATCCTCCACCAA
Reverse Sequence:	GGATTTGGGCTATGCAGGTTGG
ACCN:	<a href="#">NM_001048225</a> , <a href="#">NM_001048225.1</a> , <a href="#">NM_001048225.2</a> , <a href="#">BC001105</a> , <a href="#">BC012818</a> , <a href="#">BG703352</a> , <a href="#">BI666711</a> , <a href="#">BM767886</a> , <a href="#">BP308674</a> , <a href="#">BP330715</a> , <a href="#">BQ787725</a> , <a href="#">BX335176</a> , <a href="#">NM_001048225.4</a>
UniProt ID:	<a href="#">Q9BQY9</a>
Synonyms:	C20orf35; CK1BP; HSMNP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)