

## Product datasheet for **HP229363**

### F box protein 38 (FBXO38) Human qPCR Primer Pair (NM\_205836)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81545
Forward Sequence:	GTCCAAGAAGATGGAGAGGTGG
Reverse Sequence:	ACTACACGCTGAAGACCACTGG
ACCN:	<a href="#">NM_205836</a> , <a href="#">NM_205836.1</a> , <a href="#">BC033454</a> , <a href="#">BC033454.1</a> , <a href="#">BC005849</a> , <a href="#">BC005873</a> , <a href="#">BC015366</a> , <a href="#">BC047743</a> , <a href="#">BC050424</a> , <a href="#">BC056147</a> , <a href="#">BC078159</a> , <a href="#">BI459318</a> , <a href="#">NM_205836.3</a>
UniProt ID:	<a href="#">Q6PIJ6</a>
Synonyms:	Fbx38; HMN2D; MOKA; SP329
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)