

Product datasheet for **HP229384**

Encephalopsin (OPN3) Human qPCR Primer Pair (NM_014322)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23596
Forward Sequence:	ATGGTCACCTGGTCACTCCAAC
Reverse Sequence:	GAGGCACAGAAGCTGCAAAAGG
ACCN:	<u>NM_014322</u> , <u>NM_014322.1</u> , <u>NM_014322.2</u> , <u>BC036773</u> , <u>BC036773.1</u> , <u>BC051900</u> , <u>BC062622</u> , <u>BC073872</u> , <u>BM679966</u> , <u>NM_014322.3</u>
UniProt ID:	<u>Q9H1Y3</u>
Synonyms:	ECPN; PPP1R116
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)