

## Product datasheet for **HP229449**

### AP2M1 Human qPCR Primer Pair (NM\_004068)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1173
Forward Sequence:	GCTTTATCCCGCCAGATGGAGA
Reverse Sequence:	CTTGACCTCCAGTTTGGTGCGT
ACCN:	<a href="#">BC013796</a> , <a href="#">BC014030</a> , <a href="#">NM_004068</a> , <a href="#">NM_004068.1</a> , <a href="#">NM_004068.2</a> , <a href="#">NM_004068.3</a> , <a href="#">BC013796.1</a> , <a href="#">BC014030.1</a> , <a href="#">BC004996</a> , <a href="#">BI461415</a> , <a href="#">BT007308</a> , <a href="#">NM_004068.4</a>
UniProt ID:	<a href="#">Q96CW1</a>
Synonyms:	AP50; CLAPM1; MRD60; mu2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)