

Product datasheet for **HP229472**

FSD1L Human qPCR Primer Pair (NM_001145313)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	83856
Forward Sequence:	AAGGATGAGCGATGCTGGGAGA
Reverse Sequence:	TACTCTCCAGCCACAGCCTTGT
ACCN:	<u>BC036746</u> , <u>NM_001145313</u> , <u>NM_001145313.1</u> , <u>NM_001145313.2</u> , <u>BP219536</u> , <u>BP229804</u> , <u>BP365251</u> , <u>NM_001145313.3</u>
UniProt ID:	<u>Q9BXM9</u>
Synonyms:	CCDC10; CSDUFD1; FSD1CL; FSD1NL; MIR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)