

Product datasheet for **HP229495**

C10orf68 (CCDC7) Human qPCR Primer Pair (NM_145023)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79741
Forward Sequence:	CGGTATGCTTTGCCATTCCATC
Reverse Sequence:	CGCAATGTCCATAGGTTTCTTCC
ACCN:	<u>NM_145023</u> , <u>NM_145023.4</u> , <u>NM_145023.5</u> , <u>BC022020</u> , <u>BC042672</u> , <u>BC125090</u> , <u>BC125091</u> , <u>BI462512</u> , <u>BX647771</u> , <u>NM_145023.6</u>
UniProt ID:	<u>Q96M83</u>
Synonyms:	BIOT2; BioT2-A; BioT2-B; BioT2-C; C10orf68
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)