

## Product datasheet for **HP229514**

### **GOLPH2 (GOLM1) Human qPCR Primer Pair (NM\_177937)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	51280
<b>Forward Sequence:</b>	ATCACCACAGGTGAGAGGCTCA
<b>Reverse Sequence:</b>	ACTTCCTCTCCAGTTGGTCTG
<b>ACCN:</b>	<a href="#">NM_177937</a> , <a href="#">NM_177937.1</a> , <a href="#">NM_177937.2</a> , <a href="#">BC001740</a> , <a href="#">BC001740.1</a> , <a href="#">BQ216766</a>
<b>UniProt ID:</b>	<a href="#">Q8NBJ4</a>
<b>Synonyms:</b>	bA379P1.3; C9orf155; GOLPH2; GP73; HEL46; PSEC0257
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)