

## Product datasheet for **HP229551**

### CSPP1 Human qPCR Primer Pair (NM\_024790)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79848
Forward Sequence:	AGGCGTCTACAAGAGCGATTGC
Reverse Sequence:	CTGTCTTCTGACAGGAGCTTGC
ACCN:	<a href="#">NM_024790</a> , <a href="#">NM_024790.1</a> , <a href="#">NM_024790.2</a> , <a href="#">NM_024790.3</a> , <a href="#">NM_024790.4</a> , <a href="#">NM_024790.5</a> , <a href="#">NM_024790.6</a> , <a href="#">BC022867</a> , <a href="#">BC029445</a> , <a href="#">BC148417</a> , <a href="#">BC152993</a> , <a href="#">BX649147</a>
UniProt ID:	<a href="#">Q1MSJ5</a>
Synonyms:	CSPP; CSPP-L; JBTS21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)