

Product datasheet for **HP229552**

TMEM176B Human qPCR Primer Pair (NM_001101312)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28959
Forward Sequence:	GCGAAGTCAAGAGAACCAATGGC
Reverse Sequence:	CGAAGACCTACTCCCAAGGAAAC
ACCN:	<u>NM_001101312</u> , <u>NM_001101312.1</u> , <u>BC014586</u> , <u>BC004358</u> , <u>BC091471</u> , <u>BF511290</u> , <u>BI759037</u> , <u>BI761752</u> , <u>BI910761</u> , <u>BP298052</u> , <u>NM_001101312.2</u>
UniProt ID:	<u>Q3YBM2</u>
Synonyms:	LR8; MS4B2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)