

## Product datasheet for **HP229574**

### ZP3 Human qPCR Primer Pair (NM\_007155)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7784
Forward Sequence:	TGAGGTTGGACTCCACGAGTGT
Reverse Sequence:	GAATCTCTGCGCGTTAGTCCT
ACCN:	<a href="#">NM_007155</a> , <a href="#">NM_007155.1</a> , <a href="#">NM_007155.2</a> , <a href="#">NM_007155.3</a> , <a href="#">NM_007155.4</a> , <a href="#">NM_007155.5</a> , <a href="#">BC113949</a> , <a href="#">BC113949.1</a> , <a href="#">BC140332</a> , <a href="#">BC146482</a> , <a href="#">BI828655</a> , <a href="#">NM_007155.6</a>
UniProt ID:	<a href="#">P21754</a>
Synonyms:	OOMD3; Zp-3; ZP3A; ZP3B; ZPC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)