

Product datasheet for **HP229591**

CD64 (FCGR1B) Human qPCR Primer Pair (NM_001017986)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2210
Forward Sequence:	ATACAGGTGCCAGAGAGGTCTC
Reverse Sequence:	CCAGCTTATCCTTCCACGCATG
ACCN:	<u>NM_001017986</u> , <u>NM_001017986.1</u> , <u>NM_001017986.2</u> , <u>BC110416</u> , <u>BC156864</u> , <u>BM992325</u> , <u>BQ030076</u> , <u>BU682943</u>
Synonyms:	CD64b; FCG1; FcgammaRIa; FCGR1; FCGR1A; FcRI; IGFR1; IGFRB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)