

Product datasheet for **HP229597**

SORBS1 Human qPCR Primer Pair (NM_015385)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10580
Forward Sequence:	TATCAGCCTGGCAAGTCTTCCG
Reverse Sequence:	CCGTCTGATTCCTCTTCACT
ACCN:	<u>NM_015385</u> , <u>NM_015385.1</u> , <u>NM_015385.2</u> , <u>NM_015385.3</u> , <u>BC042612</u> , <u>BC152463</u> , <u>BC160134</u> , <u>BX640814</u>
UniProt ID:	<u>Q9BX66</u>
Synonyms:	CAP; FLAF2; R85FL; SH3D5; SH3P12; SORB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)