

Product datasheet for **HP229692**

BRD2 Human qPCR Primer Pair (NM_005104)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6046
Forward Sequence:	CGGCTTATGTTCTCCAAGTCTA
Reverse Sequence:	GGCAGTAGAGACTGGTAAAGGC
ACCN:	<u>NM_005104</u> , <u>NM_005104.1</u> , <u>NM_005104.2</u> , <u>NM_005104.3</u> , <u>BC000477</u> , <u>BC001885</u> , <u>BC007715</u> , <u>BC063840</u> , <u>BP271262</u> , <u>BX647233</u> , <u>BX648109</u>
UniProt ID:	<u>P25440</u>
Synonyms:	BRD2-IT1; D6S113E; FSH; FSRG1; NAT; O27.1.1; RING3; RNF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)