

Product datasheet for **HP229717**

RBM1A1 Human qPCR Primer Pair (NM_005058)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5940
Forward Sequence:	GGGAGTATGCTCCACCATCTAG
Reverse Sequence:	GGTGGAGCATAATCCCTGGTTTC
ACCN:	NM_005058 , NM_005058.1 , NM_005058.2 , NM_005058.3 , BC047768 , BC047768.1 , BC070298 , BC146853 , BU568309 , BX119894 , NM_005058.4
UniProt ID:	P0DJD3
Synonyms:	RBM; RBM1; RBM2; RBMY; RBMY1C; YRRM1; YRRM2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)