

## Product datasheet for **HP229816**

### PIGQ Human qPCR Primer Pair (NM\_148920)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9091
Forward Sequence:	GACTCTGCTCTTCACCATCCTG
Reverse Sequence:	TGATGAGGTCCACGAGCAGATG
ACCN:	<u><a href="#">BC010094</a></u> , <u><a href="#">NM_148920</a></u> , <u><a href="#">NM_148920.1</a></u> , <u><a href="#">NM_148920.2</a></u> , <u><a href="#">BC006377</a></u> , <u><a href="#">BC006377.2</a></u> , <u><a href="#">NM_148920.3</a></u>
UniProt ID:	<u><a href="#">Q9BRB3</a></u>
Synonyms:	c407A10.1; DEE77; EIEE77; GPI1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)