

## Product datasheet for **HP229819**

### Hexokinase 1 (HK1) Human qPCR Primer Pair (NM\_033497)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3098
Forward Sequence:	CTGCTGGTGAAAATCCGTAGTGG
Reverse Sequence:	GTCCAAGAAGTCAGAGATGCAGG
ACCN:	<a href="#">NM_033497</a> , <a href="#">NM_033497.1</a> , <a href="#">NM_033497.2</a> , <a href="#">BC008730</a> , <a href="#">BC008730.2</a> , <a href="#">BX406321</a>
UniProt ID:	<a href="#">P19367</a>
Synonyms:	hexokinase; HK; HK1-ta; HK1-tb; HK1-tc; HKD; HKI; HMSNR; HXK1; NEDVIBA; RP79
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)