

## Product datasheet for **HP229832**

### PP4R4 (PPP4R4) Human qPCR Primer Pair (NM\_058237)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57718
Forward Sequence:	GCTACTTCAGAAATATGCCTGCC
Reverse Sequence:	GTTCGTGAAGCCGCCTTTTGGGA
ACCN:	<a href="#">NM_058237</a> , <a href="#">NM_058237.1</a> , <a href="#">BC068491</a> , <a href="#">BC068491.1</a> , <a href="#">BC002650</a> , <a href="#">BI546744</a> , <a href="#">NM_058237.2</a>
UniProt ID:	<a href="#">Q6NUP7</a>
Synonyms:	CFAP14; KIAA1622; PP4R4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)