

Product datasheet for **HP229837**

RFC2 Human qPCR Primer Pair (NM_181471)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5982
Forward Sequence:	GTCGGGAATGAAGACACCGTGA
Reverse Sequence:	CAGAATGCTTGTGGTCTTGCCG
ACCN:	<u>NM_181471</u> , <u>NM_181471.1</u> , <u>NM_181471.2</u> , <u>BC002813</u> , <u>BC002813.1</u> , <u>BM563468</u> , <u>BM810551</u> , <u>BQ070526</u> , <u>BQ674131</u> , <u>BT007058</u> , <u>NM_181471.3</u>
UniProt ID:	<u>P35250</u>
Synonyms:	RFC40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)