

Product datasheet for **HP229846**

CYB561 Human qPCR Primer Pair (NM_001017917)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 1534 |
| Forward Sequence: | CCTGCTGGTTTACCGTGTCTTC |
| Reverse Sequence: | CTTCCTGTGGTAGTCGAACACC |
| ACCN: | NM_001017917 , NM_001017917.1 , BC091485 , BC091485.1 , BC000021 , BC002976 , BF206825 , BM552101 , BT007096 , NM_001017917.2 |
| UniProt ID: | P49447 |
| Synonyms: | CYB561A1; FRRS2; ORTHYP2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)