

Product datasheet for **HP229866**

PLSCR4 Human qPCR Primer Pair (NM_001128304)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57088
Forward Sequence:	CCTTCAGATGCACCTGCTGTTG
Reverse Sequence:	CCGCAACAAAGCCAATGGTGAC
ACCN:	<u>NM_001128304</u> , <u>NM_001128304.1</u> , <u>BC028354</u> , <u>BC028354.2</u> , <u>BG538024</u> , <u>BI598345</u> , <u>BP201782</u> , <u>BU160106</u> , <u>BU567842</u> , <u>NM_001128304.2</u>
UniProt ID:	<u>Q9NRQ2</u>
Synonyms:	TRA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)