

Product datasheet for **HP229871**

CD82 Human qPCR Primer Pair (NM_002231)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3732
Forward Sequence:	GGAAGAGGACAACAGCCTTTCTG
Reverse Sequence:	ATGATGCCCAGTTCTCCTGCA
ACCN:	<u>NM_002231</u> , <u>NM_002231.1</u> , <u>NM_002231.2</u> , <u>NM_002231.3</u> , <u>BC000726</u> , <u>BC000726.2</u> , <u>BC001821</u> , <u>BT019761</u> , <u>BT019762</u> , <u>NM_002231.4</u>
UniProt ID:	<u>P27701</u>
Synonyms:	4F9; C33; GR15; IA4; KAI1; R2; SAR2; ST6; TSPAN27
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)