

Product datasheet for **HP229876**

SLC6A6 Human qPCR Primer Pair (NM_001134367)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6533
Forward Sequence:	TGATGTGGCTGAGTCAGGTCCT
Reverse Sequence:	CTGGCTATCCAGTCCAAGCAAG
ACCN:	<u>NM_001134367</u> , <u>NM_001134367.1</u> , <u>NM_001134367.2</u> , <u>NM_001134367.3</u> , <u>BC006252</u> , <u>BC111489</u> , <u>BC137128</u> , <u>BC137129</u> , <u>BC144490</u> , <u>BM149397</u> , <u>BM873716</u> , <u>BM981004</u> , <u>BQ945675</u> , <u>BU625571</u> , <u>BX478649</u> , <u>BX509263</u> , <u>BX510251</u>
UniProt ID:	<u>P31641</u>
Synonyms:	TAUT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)