

Product datasheet for **HP229878**

SLC5A9 Human qPCR Primer Pair (NM_001135181)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	200010
Forward Sequence:	TTCCGCAGGAAGTCAACAGAGC
Reverse Sequence:	GGTAACTGGTGACAGCCTGGAT
ACCN:	<u>NM_001135181</u> , <u>NM_001135181.1</u> , <u>BC025411</u> , <u>BC032417</u> , <u>BC048332</u> , <u>BC104857</u> , <u>BC104863</u> , <u>BC143570</u> , <u>BC143572</u> , <u>BC143587</u> , <u>BX648549</u> , <u>NM_001135181.2</u>
UniProt ID:	<u>Q2M3M2</u>
Synonyms:	SGLT4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)