

## Product datasheet for **HP229886**

### GPR126 (ADGRG6) Human qPCR Primer Pair (NM\_198569)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57211
Forward Sequence:	ACAGAGCAAGGTGGCAGAATGG
Reverse Sequence:	TTGTCCTCTCCAGCACTCAGGT
ACCN:	<u><a href="#">NM_198569</a></u> , <u><a href="#">NM_198569.1</a></u> , <u><a href="#">NM_198569.2</a></u> , <u><a href="#">BC075798</a></u> , <u><a href="#">BC075798.1</a></u> , <u><a href="#">BC036008</a></u> , <u><a href="#">BX640873</a></u> , <u><a href="#">BX640971</a></u> , <u><a href="#">BX648315</a></u>
UniProt ID:	<u><a href="#">Q86SQ4</a></u>
Synonyms:	APG1; DREG; GPR126; LCCS9; PR126; PS1TP2; STQTL1; VIGR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)