

Product datasheet for **HP229892**

SLC22A23 Human qPCR Primer Pair (NM_021945)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63027
Forward Sequence:	CCCAAGAAGGTCTGCATCGTGA
Reverse Sequence:	TCCTGGCAAAGCAGTGGTGGAT
ACCN:	<u>NM_021945</u> , <u>NM_021945.3</u> , <u>NM_021945.4</u> , <u>NM_021945.5</u> , <u>NM_021945.6</u> , <u>BC128580</u> , <u>BC022217</u> , <u>BC038748</u> , <u>BC040876</u> , <u>BC090851</u> , <u>BC128581</u>
UniProt ID:	<u>A1A5C7</u>
Synonyms:	C6orf85
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)