

Product datasheet for **HP229917**

SPOCK3 Human qPCR Primer Pair (NM_016950)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	50859
Forward Sequence:	GTGGTCTATCCCAGCCCTGTTT
Reverse Sequence:	GGGCTTATCTGAAGGACATGGG
ACCN:	<u>NM_016950</u> , <u>NM_016950.1</u> , <u>NM_016950.2</u> , <u>BC013983</u> , <u>BC013983.2</u> , <u>BC000460</u> , <u>BC003017</u> , <u>BI603782</u>
UniProt ID:	<u>Q9BQ16</u>
Synonyms:	HSAJ1454; TES-3; TICN3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)