

Product datasheet for **HP229924**

GTPBP9 (OLA1) Human qPCR Primer Pair (NM_013341)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29789
Forward Sequence:	GACAGAAGTATCTGGAAGCGAAC
Reverse Sequence:	CTGATGGTCCATGCACGCACTT
ACCN:	<u>NM_013341</u> , <u>NM_013341.1</u> , <u>NM_013341.2</u> , <u>NM_013341.3</u> , <u>NM_013341.4</u> , <u>BC012842</u> , <u>BC013925</u> , <u>BC029376</u> , <u>BC068447</u> , <u>BC091522</u> , <u>NM_013341.5</u>
UniProt ID:	<u>Q9NTK5</u>
Synonyms:	DOC45; GBP45; GTBP9; GTPBP9; PTD004
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)