

Product datasheet for **HP229937**

DLX4 Human qPCR Primer Pair (NM_138281)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1748
Forward Sequence:	GCAGCACCTAAACCAGCGTTTC
Reverse Sequence:	ACTTGAGCGTTTGTCTGAAACC
ACCN:	<u>NM_138281</u> , <u>NM_138281.1</u> , <u>NM_138281.2</u> , <u>BC016145</u> , <u>BC016145.1</u> , <u>BC005812</u> , <u>BC014419</u> , <u>BG479840</u> , <u>BT006978</u> , <u>BX105774</u> , <u>BX113580</u> , <u>NM_138281.3</u>
UniProt ID:	<u>Q92988</u>
Synonyms:	BP1; DLX7; DLX8; DLX9; OFC15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)