

Product datasheet for **HP229992**

RABEP1 Human qPCR Primer Pair (NM_004703)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9135
Forward Sequence:	AGCAGCGAAGATTCGAGTCACC
Reverse Sequence:	CTGAACATCCCTCTTGCCTGG
ACCN:	NM_004703 , NM_004703.1 , NM_004703.2 , NM_004703.3 , NM_004703.4 , NM_004703.5 , BC041700 , BC041700.1 , BC018742 , BC032435 , BC065242 , BI666786 , BM678290 , BX103338
UniProt ID:	Q15276
Synonyms:	RAB5EP; RABPT5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)