

Product datasheet for **HP230013**

FAM178A (SLF2) Human qPCR Primer Pair (NM_001136123)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55719
Forward Sequence:	CCTACTTCAGCAGGCTCTAAGC
Reverse Sequence:	TGGTCCATCAGGACGCTTCTCA
ACCN:	BC015131 , NM_001136123 , NM_001136123.1 , BC020915 , BC030565 , BC047715 , BC064526 , BC073832 , BC160007 , NM_001136123.2
UniProt ID:	Q8IX21
Synonyms:	C10orf6; FAM178A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)