

Product datasheet for **HP230060**

Artemis (DCLRE1C) Human qPCR Primer Pair (NM_022487)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 64421 |
| Forward Sequence: | GAGCTAGAACAGTTCACCGAGAC |
| Reverse Sequence: | CAGGCTGCTTTTCTGATACTGCA |
| ACCN: | NM_022487 , NM_022487.1 , NM_022487.2 , NM_022487.3 , BC009185 , BC009185.2 , BC000863 , BC022254 , BC067233 , BC108276 , BE799340 , BQ932042 , BX648858 , NM_022487.4 |
| UniProt ID: | Q96SD1 |
| Synonyms: | A-SCID; DCLREC1C; RS-SCID; SCIDA; SNM1C |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)