

Product datasheet for **HP230061**

DPH2 Human qPCR Primer Pair (NM_001384)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1802
Forward Sequence:	GATCATGGAAGCTTGGCTCTGAC
Reverse Sequence:	ACAGCTTCTGTCACTGGCGTCT
ACCN:	NM_001384 , NM_001384.1 , NM_001384.2 , NM_001384.3 , NM_001384.4 , BC003181 , BC003181.1 , BC001389 , BC016956 , BG385934 , BG472282 , BG747150 , BM564151 , BM794530 , BQ056213 , BQ215136 , BT007431 , BX362196 , NM_001384.5
UniProt ID:	Q9BQC3
Synonyms:	DPH2L2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)