

Product datasheet for **HP230089**

Calcium independent Phospholipase A2 (PLA2G6) Human qPCR Primer Pair (NM_003560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8398
Forward Sequence:	CCTCATCGTGTTCCGAGCAGAA
Reverse Sequence:	CACGGTTCTCAGCAGAGTCAAG
ACCN:	<u>NM_003560</u> , <u>NM_003560.1</u> , <u>NM_003560.2</u> , <u>NM_003560.3</u> , <u>BC036742</u> , <u>BC036742.1</u> , <u>BC051904</u> , <u>BC051904.1</u> , <u>BC034592</u> , <u>NM_003560.4</u>
UniProt ID:	<u>O60733</u>
Synonyms:	CaI-PLA2; GVI; INAD1; iPLA2; IPLA2-VIA; iPLA2beta; NBIA2; NBIA2A; NBIA2B; PARK14; PLA2; PNPLA9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)