

Product datasheet for **HP230092**

HEY1 Human qPCR Primer Pair (NM_012258)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23462
Forward Sequence:	TGTCTGAGCTGAGAAGGCTGGT
Reverse Sequence:	TTCAGGTGATCCACGGTCATCTG
ACCN:	<u>NM_012258</u> , <u>NM_012258.1</u> , <u>NM_012258.2</u> , <u>NM_012258.3</u> , <u>BC001873</u> , <u>BC001873.1</u> , <u>BP283540</u> , <u>BT020065</u> , <u>NM_012258.4</u>
UniProt ID:	<u>Q9Y5J3</u>
Synonyms:	BHLHb31; CHF2; HERP2; HESR1; hHRT1; HRT-1; NERP2; OAF1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)