

Product datasheet for **HP230099**

Growth Arrest Specific Protein 7 (GAS7) Human qPCR Primer Pair (NM_201433)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8522
Forward Sequence:	CTCTCAGAACTCCTTGGCTTCAC
Reverse Sequence:	GTTCTCACGGAAGTTCATCAGGG
ACCN:	<u>NM_201433</u> , <u>NM_201433.1</u> , <u>BC006454</u> , <u>BC006454.1</u> , <u>BC001152</u> , <u>BC113835</u> , <u>BP202727</u> , <u>BT006891</u> , <u>NM_201433.2</u>
UniProt ID:	<u>O60861</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)