

Product datasheet for **HP230117**

CLDND1 Human qPCR Primer Pair (NM_001040182)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56650
Forward Sequence:	TGGTATAGCCCACCAGAAAGGAC
Reverse Sequence:	CAATCCCGCTATTGTGGTTCCG
ACCN:	<u>NM_001040182</u> , <u>NM_001040182.1</u> , <u>BC001757</u> , <u>BC001757.2</u> , <u>BC022551</u> , <u>BC022551.1</u> , <u>BC013610</u> , <u>BC095441</u> , <u>BF967477</u> , <u>BI552824</u> , <u>BU784186</u>
UniProt ID:	<u>Q9NY35</u>
Synonyms:	C3orf4; GENX-3745; Z38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)