

Product datasheet for **HP230168**

SERPINB7 Human qPCR Primer Pair (NM_003784)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8710
Forward Sequence:	CGATGCCAAAGTGGAGCGAGTT
Reverse Sequence:	CAGCAGATGAGCTTATGCCACC
ACCN:	NM_003784 , NM_003784.1 , NM_003784.2 , NM_003784.3 , BC069417 , BC069442 , BC069547 , BC106743 , BC106744 , BP332064 , BX641156 , NM_003784.4
UniProt ID:	O75635
Synonyms:	MEGSIN; PPKN; TP55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)