

Product datasheet for **HP230221**

GLEPP1 (PTPRO) Human qPCR Primer Pair (NM_030667)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5800
Forward Sequence:	CTGCTGTGACAGCTCTACCATC
Reverse Sequence:	GCAACCACCATCCAGTGAAGCA
ACCN:	<u>NM_030667</u> , <u>NM_030667.1</u> , <u>NM_030667.2</u> , <u>BC126201</u> , <u>BC035960</u> , <u>BC126203</u> , <u>NM_030667.3</u>
UniProt ID:	<u>Q16827</u>
Synonyms:	GLEPP1; NPHS6; PTP-OC; PTP-U2; PTPROT; PTPU2; R-PTP-O
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)