

Product datasheet for **HP230236**

SLC25A45 Human qPCR Primer Pair (NM_182556)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283130
Forward Sequence:	ACCCGTTTGACACTGTAAAGGTG
Reverse Sequence:	CAGGACAGAGTTGACCACAGCT
ACCN:	<u>BC036869</u> , <u>NM_182556</u> , <u>NM_182556.1</u> , <u>NM_182556.2</u> , <u>NM_182556.3</u> , <u>BC041100</u> , <u>BM695511</u> , <u>BU728924</u> , <u>BX648573</u>
UniProt ID:	<u>Q8N413</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)