

Product datasheet for **HP230267**

PAICS Human qPCR Primer Pair (NM_006452)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10606
Forward Sequence:	ACCACCTGGAAGGAAAAGCTGC
Reverse Sequence:	CGGTGCAATGAAAGCTGTCTCC
ACCN:	NM_006452 , NM_006452.1 , NM_006452.2 , NM_006452.3 , BC010273 , BC019255 , BM807012 , BT006988
UniProt ID:	P22234
Synonyms:	ADE2; ADE2H1; AIRC; PAIS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)