

Product datasheet for **HP230279**

BCOR Human qPCR Primer Pair (NM_017745)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54880
Forward Sequence:	TGTCTACCCGCTGCTTACTGTG
Reverse Sequence:	TCTCGGAGTCTTTGGTTGCTGG
ACCN:	NM_017745 , NM_017745.1 , NM_017745.2 , NM_017745.3 , NM_017745.4 , NM_017745.5 , BC114220 , BC009675 , BC063536 , BC128387 , BC128456
UniProt ID:	Q6W2I9
Synonyms:	ANOP2; MAA2; MCOPS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)