

Product datasheet for **HP230303**

GGPS1 Human qPCR Primer Pair (NM_004837)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9453
Forward Sequence:	TGAGAACCAAACCTTTCACAGGCAT
Reverse Sequence:	CGATGAGTAAACTGGCATTATGCA
ACCN:	BC032569 , NM_004837 , NM_004837.2 , NM_004837.3 , BC005252 , BC005252.1 , BC067768 , BI462253 , BI520491 , BI963107 , NM_004837.4
UniProt ID:	O95749
Synonyms:	GGPPS; GGPPS1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)