

Product datasheet for **HP230311**

FRMPD2 Human qPCR Primer Pair (NM_152428)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	143162
Forward Sequence:	ACTTCAACGTGGACGGCAGCAA
Reverse Sequence:	GGTCATTGGCTTCTGGATCACAC
ACCN:	<u>NM_152428</u> , <u>NM_152428.4</u> , <u>NM_152428.1</u> , <u>BC031614</u> , <u>BC040910</u> , <u>BC073954</u> , <u>BC144636</u> , <u>BC144639</u> , <u>BC156408</u> , <u>BC157131</u>
UniProt ID:	<u>Q68DX3</u>
Synonyms:	FERM and PDZ domain containing 2; MGC35285; MGC87776; MGC87777; MGC90186; PDZD5C; PDZ domain containing 5C; PDZK5C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)